HYALUAL®

INSTRUCTIONS FOR USE

Injectable implant based on hyaluronic acid Hyalual®**Description**

Hyalual® is an elastic gel based on hyaluronic acid of non-animal origin. A sterile, non-pyrogenic and transparent gel.

Contents of package

There is 1 pre-filled syringe with a **Hyalual®** implant in each package. There are two sterile 30G needles for each syringe in the package. The volume of each syringe is shown on both the cardboard package and on each syringe.

Composition

1ml volume of **Injectable implant based on hyaluronic acid Hyalual® light (1.1%)** contains:

Sodium Hyaluronate.....11.0mg;
Succinate buffer pH 7.4.....up to 1.0ml.

2ml volume of **Injectable implant based on hyaluronic acid Hyalual® light (1.1%)** contains:

Sodium Hyaluronate.....22.0mg;
Succinate buffer pH 7.4.....up to 2.0ml.

1ml volume of **Injectable implant based on hyaluronic acid Hyalual® smart (1.8%)** contains:

Sodium Hyaluronate.....18.0mg;
Succinate buffer pH 7.4.....up to 1.0ml.

2ml volume of **Injectable implant based on hyaluronic acid Hyalual® smart (1.8%)** contains:

Sodium Hyaluronate.....36.0mg;
Succinate buffer pH 7.4.....up to 2.0ml.

1ml volume of **Injectable implant based on hyaluronic acid Hyalual® booster (2.2%)** contains:

Sodium Hyaluronate.....22.0mg;
Succinate buffer pH 7.4.....up to 1.0 ml.

2ml volume of **Injectable implant based on hyaluronic acid Hyalual® booster (2.2%)** contains:

Sodium Hyaluronate.....44.0mg;
Succinate buffer pH 7.4.....up to 2.0ml.

Indications for use

- for use in reconstructive surgery by filling deformations caused by injuries, and restoring volumes lost as a result of age changes;
- changes the skin anatomy by restoring lost volumes and skin hydration, filling folds and hollows of skin.

Action mechanism

Hyalual® is a biologically compatible long-action injection implant. Thanks to its especially viscoelastic properties, it creates a natural, uniform and homogeneous volume when adjusting skin that has lost its tone or become dehydrated or when filling folds and hollows of skin. The hyaluronic acid contained in **Hyalual®** guarantees long-term absorption of the injection implant by the skin tissues. This characteristic allows support of the optimal level of derma hydration, assisting the prolonged elasticity of skin. In order to obtain optimal results, the need to repeat one or several sessions may arise, depending on the recommendations of the specialist conducting the procedures.

Contraindications

Hyalual® should not be used:

- for other forms of injection besides intradermal;
- in combination with peeling, laser or ultrasound procedures;
- if the patient has skin diseases, infections, or infected areas in the treatment area or around it;
- if the patient is known to be hypersensitive to hyaluronic acid, or has a history of allergies or anaphylactic shock;
- in the case of the patient having autoimmune diseases;
- during pregnancy and lactation;
- for children;
- to inject into blood vessels;
- to inject into the eyelid area.

As interaction with other well-known injection products has not been researched, it is not recommended to use **Hyalual®** jointly with other injection products around or in the treatment area.

Method of use and doses

- Injection of a **Hyalual®** implant must be carried out by a specialist in possession of the technology to inject injection implants.
- Before the start of procedures, detailed information must be obtained from the patient about his/her former illnesses and the patient informed about the possible results of the implant and possible side effects.
- Carefully disinfect the treatment area.
- The use of the needles contained in the package is recommended for injection of the implant.
- If necessary, local and regional anaesthetic may be used. In such a case, the instructions for use for these treatments must be followed.
- Make sure that the needle is firmly connected to the syringe. If there is a blockage in the needle or complication in injecting the implant, halt the procedure and change the needle.
- The choice of technique and the number of procedures depends on the type of the patient's skin, the characteristics of the defect being corrected, and the method of injection.

Hyalual® is injected into the middle layer of the dermis in the dehydrated area and/or the folds and hollows requiring correction. If the implant is injected too deep, i.e. into a deep layer of the dermis or the subcutaneous tissues, the correction will be less effective.

- The volume of the implant injected depends on the degree of correction required. It is important not to make overcorrection.
- As a rule, for optimal elimination of visible aesthetic defects, 3-5 procedures with an interval of 1-2 weeks are required. Carrying out additional procedures ensures long-term effects.

Side effects

- allergic and skin reactions: reddening, swelling, oedema, hematomas, itching, and mild pain in the injection area can appear after injection, however, as a rule, they disappear within 72 hours;
- hardening, appearance of blotches (specific colouring or discolouring of the skin in the injection area), loss of feeling in the area of injection of the implant;
- infections; cases of necrosis, formation of tubercles, granules, heightened sensitivity and purulent inflammations after an injection of hyaluronic acid have been recorded.

The specialist responsible for injection must inform a patient about the potential risks beforehand.

The patient must inform the doctor about all other side effects, as well as if the side effects described last for more than one week. The doctor in turn is obliged to inform the manufacturer or distributor as soon as possible.

Incompatibility

Injectable implant based on hyaluronic acid Hyalual® must not be used with quaternary ammonium salts such as benzalkonium chloride.

Safety measures

- Patients must not take large doses of Vitamin E, acetylsalicylic acid, antipyretics or anticoagulants for one week before the procedures are carried out.
- Particular care is recommended when injecting the product in the anatomical nerve area. Violation of the implant injection technique in these areas could cause injury to the peripheral nerves, damage sensitivity and lead to paralysis of the affected innervation areas.
- If there is a history of herpes rash, prophylactic treatment of herpes is recommended before undergoing the procedure.
- Patients are not recommended to use cosmetics for 12 hours after the injection.
- Avoid sharp temperature fluctuations (intense cold, saunas, steam baths) for one week after the procedure.
- In order to prevent the movement of the implant, patients are not recommended to touch the place of injection, and should avoid injury to the area after the injection.
- Check that the package is sealed before use.
- Check the use-by date on the label.
- Do not use an implant if the use-by date has passed, the hermetic sealing of the packet has been damaged, or the contents of the syringe have changed colour or become opaque.
- Do not re-use.
- Do not mix with other products.
- Must not be resterilised.
- Dispose of used needles and syringes in accordance with the local requirements of waste disposal.

- Do not try to straighten out a bent needle or cannula; it must be replaced.


Use-by date: 2 years.

Storage conditions

Store in a dry place away from sunlight inaccessible to children, at a temperature of between +5 °C and +30 °C. Do not freeze.

Name and address of manufacturer

Yuria-Pharm Ltd., 10, M. Amosova Str., Kyiv, Ukraine, 03038
 Tel.: +38 (044) 275-92-42; +38 (044) 275-01-08 www.uf.ua
 Production facility address:
 108, Kobzarska Str., Cherkasy, Ukraine, 18030

Graphical symbols and their interpretation	
	STERILIZED USING STEAM OR DRY HEAT
	STERILE MEDICAL DEVICE IN PRIMARY PACKAGING
	DO NOT RE-USE
	MANUFACTURER
	KEEP AWAY FROM SUNLIGHT
	KEEP DRY
	DO NOT RESTERILIZE
	DO NOT USE IF PACKAGE IS DAMAGED
	FRAGILE, HANDLE WITH CARE
	TEMPERATURE LIMIT
	CAUTION
	CONSULT INSTRUCTIONS FOR USE
	DATE OF MANUFACTURE
	USE-BY DATE
	BATCH CODE
	NON-PYROGENIC
	MARK OF COMPLIANCE WITH DIRECTIVE 93/42/EEC ON MEDICAL DEVICES AND CODE OF THE BODY VERIFYING COMPLIANCE
	CONFORMITY MARK TO THE TECHNICAL REGULATION FOR MEDICAL DEVICES AND THE CODE OF THE CONFORMITY ASSESSMENT BODY
	AUTHORIZED REPRESENTATIVE IN THE EUROPEAN COMMUNITY

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